

Data sheet

BL 250 kit 015 and 025 CVVH/CVVHD



Product description

The procedure packs are captive to use on the CARPEDIEM™ machine. These procedure packs are composed by BL 250 preassembled device for hemofiltration/hemodialysis 015 for CARPEDIEM™ machine or BL 250 preassembled device for hemofiltration/hemodialysis 025 for CARPEDIEM™ machine manufactured by Bellco S.r.l and a 10 ml Luer Lok™ BD Plastipak syringe, that is manufactured by Becton Dickinson S.A.

Each preassembled device is sterile, single-use, extracorporeal circuit and consists of a dialyzer and tubing lines, that are permanently connected, and three 3-liter waste bags. The dialyzer consists of a cylindrical body that contains a bundle of hollow fibers made of high permeability polyethersulfone.

The 10 ml Luer Lok™ BD Plastipak syringe is sterile, single-use and inserted with the primary packaging in each procedure pack.

The two preassembled devices have different hollow-fiber dialyzer dimensions, HCD 015 and HCD 025, but they have the same tubing line configuration: a blood access line, a blood return line, a dialysate/infusion line, an effluent outlet line, and a heparin infusion line (figure 1).




Figure 1 - legend

- | | | |
|-----------------------------|-------------------------------------|--------------------------------------|
| 1. Arterial line Luer-lock | 7. Venous chamber | 13. Heparin pump line Luer-lock |
| 2. Arterial pressure dome | 8. Venous line Luer-lock | 14. Venous return line |
| 3. Blood pump segment | 9. Infusion/dialysis line Luer-lock | 15. Effluent line |
| 4. Pre-filter pressure dome | 10. Infusion/dialysis pump segment | 16. Arterial line sampling port |
| 5. Filter | 11. Effluent pump segment | 17. Venous return line sampling port |
| 6. Venous pressure dome | 12. Effluent line Luer-lock | 18. Heparin line |

Applicable therapies

The procedure packs for CARPEDIEM™ machine consist of a preassembled device and a 10 ml Luer Lok™ BD Plastipak syringe. The therapies described below depend on the use of these preassembled devices that form the procedure packs with the CARPEDIEM™ machine.

The preassembled devices are indicated for use with the CARPEDIEM™ machine in extracorporeal blood treatment for pediatric patients weighing 2.5 kilograms or more with acute kidney injury, fluid overload and/or electrolyte disorders, requiring, hemofiltration (CVVH), hemodialysis (CVVHD) or ultrafiltration/fluid removal (SCUF).

Depending on the configuration, these devices can perform the following treatments:

- Continuous Venous-Hemodialysis (CVVHD) a continuous form of hemodialysis characterized by a slow dialysate flow (in-flow or counter-flow) with respect to the blood flow in the dialysate compartment of the hemodialyzer. The main solute transmembrane removal mechanism is diffusion.
- Continuous Venous-Hemofiltration (CVVH) a continuous form of hemofiltration with infusion of replacement fluid upstream (pre-infusion) or downstream (post-infusion) from the hemofilter. The solute transmembrane transport mechanism is convection. The ultrafiltrate is partly or completely replaced with an appropriate replacement fluid in order to obtain effective reduction of the solute concentration and effective fluid balance.
- Slow Continuous Ultrafiltration (CVVH-SCUF) is a treatment based solely on slow removal of plasma water. It is generally used to manage patients who have fluid overload and/or are immune to pharmacologic treatment and may also have concomitant renal damage. The goal of the treatment is to achieve effective and safe correction of fluid overload. The solute concentration cannot be reduced using this treatment.

Intended use

Procedure packs for CVVH (Continuous Venous-Hemofiltration) and CVVHD (Continuous Venous-Hemodialysis) therapies.

Codes available

| Procedure pack REF | Procedure pack name | Medical device REF | Medical device name | Medical device manufacturer | CND ² | GMDN ² |
|-----------------------|-----------------------------|-----------------------|--|--|------------------|-------------------|
| IB0595540 | BL250 KIT 015 CVVH/CVVHD | IB0580804 | BL250 Preassembled device for hemofiltration/hemodialysis 015 for CARPEDIEM™ machine | Bellco S.r.l. - Via Camurana 1 - 41037 Mirandola (MO) Italy. CE0123 | F0306 | 61674 |
| | | 305959 | BD Plastipak | Becton Dickinson S.A., Camino de Valdeoliva, s/n 28750 San Agustín del Guadalix, Madrid, Spain. CE0318 | A020102020102 | 47017 |
| IB0595550 | BL250 KIT 025 CVVH/CVVHD | IB0580805 | BL250 Preassembled device for hemofiltration/hemodialysis 025 for CARPEDIEM™ machine | Bellco S.r.l. - Via Camurana 1 - 41037 Mirandola (MO) Italy. CE0123 | F0306 | 61674 |
| | | 305959 | BD Plastipak | Becton Dickinson S.A., Camino de Valdeoliva, s/n 28750 San Agustín del Guadalix, Madrid, Spain. CE0318 | A020102020102 | 47017 |

² The CND and GMDN indicated in the table refer to the preassembled and Luer Lok™ BD Plastipak syringe.

Sterilization method and validity

The Procedure Packs are non sterile and non-pyrogenic and their shelf life is of 3 years, consist of a preassembled device and a 10 ml Luer Lok™ BD Plastipak syringe.

The shelf life of the Procedure Packs is determined by the component with the shortest shelf life.

In particular, the preassembled device and the syringe are sterile and non-pyrogenic, their sterilizing agent is ethylene oxide and they cannot be re-sterilized. The shelf life of preassembled device is 3 years and instead of the syringe is 5 years.

Technical characteristics

The technical characteristics of the procedure packs components for CARPEDIEM™ machine are reported below.

| Dialyzer | |
|------------|------------------|
| Components | Materials |
| Membrane | Polyethersulfone |
| Housing | Copolyester |
| Header | Copolyester |
| Potting | Polyurethane |
| O-ring | Silicone rubber |

| Blood lines | |
|--------------------------------------|---|
| Components | Materials |
| Tubing lines | Polyvinyl chloride |
| Blood pump segment | Polyvinyl chloride |
| Infusion pump segment | Polyvinyl chloride |
| Heparin line | Polyvinyl chloride |
| Pressure transducer membrane | Silicon rubber |
| Tube adapter | Polyvinyl chloride |
| Line connector | Polyvinyl chloride (PVC) - Methyl methacrylate-acrylonitrile-butadiene-styrene (MABS) |
| Pressure traducer holder | Polyvinyl chloride |
| Venous chamber | Polyvinyl chloride |
| Vonous chamber filter | Polyethylene |
| Filter connector | Polyvinyl chloride |
| Access port | Polyvinyl chloride- Isoprene - Polypropylene |
| Clamps | Polypropylene |
| INF/UF luer connector ring | Polycarbonate |
| Vonous drip chamber service line cap | Polyvinyl chloride (PVC) |
| INF Y connector | Polyvinyl chloride |
| One-way valve | Silicon rubber - Methyl methacrylate-acrylonitrile-butadiene-styrene (MABS) |
| Port caps | Polypropylene |
| Vented spike | Acrylonitrile-butadiene-styrene (ABS) - Low density polyethylene (LDPE) - Linear Low density polyethylene (LLDPE) |
| Unvented spike | Acrylonitrile butadiene styrene - Low density polyethylene |
| INF/UF luer connector | Polyvinyl chloride |

| 3-liter waste bag | |
|------------------------|-----------------------------|
| Components | Materials |
| Film | Polyvinyl chloride DOP free |
| Tube | Polyvinyl chloride DOP free |
| Joint female luer lock | Polyvinyl chloride |
| Male luer lock caps | Polypropylene |
| Clamp | Polypropylene |

| Syringe ³ | |
|----------------------|--|
| Components | Materials |
| Barrel | Polypropylene |
| Piston valve | Polypropylene |
| Stoppers | Elastomer (LATEX free) |
| Lubricant | Silicone Oil (<0,25 mg/cm ²) |

³ See the technical data sheet (section n°13) provides by the manufacturer by Becton Dickinson S.A.

| Technical characteristics | | | | | | | | | | | | |
|---------------------------|--|------------------------------------|---------------------------------------|--|-----------------------------|--------------------------------------|--|--|--------------------------|---|---|---|
| Model | Dialyzer Surface area ⁴ (m ²) | Dialyzer fiber wall thickness (µm) | Dialyzer fiber internal diameter (µm) | Dialyzer blood compartment priming volume (ml) | Dialyzer maximum TMP (mmHg) | Dialyzer maximum blood flow (ml/min) | Dialyzer maximum dialysate flow (ml/min) | Dialyzer blood compartment pressure drop ⁵ (mmHg) | | Dialyzer dialysis fluid compartment pressure drop ⁶ (mmHg) | Dialyzer total length ⁷ (mm) | Total priming volume preassembled device (ml) |
| | | | | | | | | Q _b 10 ml/min | Q _b 50 ml/min | | | |
| HCD 015 | 0.17 | 30 | 200 | 11 | 500 | 50 | 10 | 19 | 32 | 10 | 128 | 32 |
| HCD 025 | 0.29 | 30 | 200 | 20 | 500 | 50 | 10 | 22 | 35 | 17 | 140 | 41 |

⁴ Medium value ±10% in according with IFU

⁵ Bovine blood: Hct= 32±3%, proteins = 60±5 g/l

⁶ Dialysis fluid: NaCl = 0.9%

⁷ Dialyzer outer body characteristic

Performance

The measurements in the charts below are taken in accordance with EN ISO 8637-1⁸. The values indicated are to be considered approximate and may vary due to measurement methods, inherent variations of the membrane, manufacturing and storage conditions. During the treatment, performance on the individual patient may vary due to variable clinical parameters of the patient.

⁸ The different ISO standard mentioned in the document are in current revision.

The performances indicated below refer to the preassembled device

| Dialyzer In Vitro clearance ⁹ | | | | | | | | | | | | | | | | |
|--|------------------|-------------------|-------------------|-------------------|---------------------|-------------------|-------------------|-------------------|---------------------|-------------------|-------------------|-------------------|----------------------|-------------------|-------------------|-------------------|
| Model | Urea (ml/min) | | | | Creatinine (ml/min) | | | | Phosphates (ml/min) | | | | Vitamin B12 (ml/min) | | | |
| | Q _b 5 | Q _b 10 | Q _b 20 | Q _b 50 | Q _b 5 | Q _b 10 | Q _b 20 | Q _b 50 | Q _b 5 | Q _b 10 | Q _b 20 | Q _b 50 | Q _b 5 | Q _b 10 | Q _b 20 | Q _b 50 |
| HCD 015 | 2.8 | 4.2 | 5.6 | 7.8 | 2.9 | 4.3 | 5.7 | 7.8 | 3.0 | 4.4 | 5.5 | 7.7 | 3.0 | 4.3 | 5.0 | 6.5 |
| HCD 025 | 2.6 | 4.0 | 5.4 | 7.2 | 2.8 | 4.2 | 5.8 | 7.5 | 2.8 | 4.3 | 5.8 | 7.2 | 2.9 | 4.3 | 5.5 | 6.7 |

⁹ In vitro clearance: Q_b = 10 ml/min, Q_f = 0 ml/min; saline solution: NaCl = 0.9%

| Clearance at maximum Q _f and Q _b (ml/min) ¹⁰ | | | | |
|---|------|------------|------------|-------------|
| Dialyzer | Urea | Creatinine | Phosphates | Vitamin B12 |
| HCD 015 ¹¹ | 10.6 | 10.7 | 10.6 | 9.2 |
| HCD 025 ¹² | 10.9 | 11.1 | 11.1 | 10.5 |

¹⁰ Saline solution: NaCl = 0.9%, Q_b = 10 ml/min

¹¹ Q_f = 14 ml/min, Q_b = 50 ml/min, Q_d = 10 ml/min

¹² Q_f = 15 ml/min, Q_b = 50 ml/min, Q_d = 10 ml/min

| Ultrafiltration coefficient | | |
|-----------------------------|--|--------------------------------|
| Dialyzer | K _{uf} ¹³ (ml/h*mmHg) | |
| | Q _B 10 ml/min | Q _B 50 ml/min |
| HCD 015 | 4.8 | 9.8 |
| HCD 025 | 9.0 | 22.1 |

¹³ Bovine blood Hct = 25±2%, proteins = 60 g/l

| Sieving coefficient | | |
|---------------------|---------|---------|
| Marker | HCD 015 | HCD 025 |
| Inulin | 0.8 | 0.8 |
| Myoglobin | 0.34 | 0.34 |
| Albumin | 0.002 | 0.002 |

Packaging

| Model | Primary packaging | | |
|---|--|---------------|------------------|
| | Pouch material | Tray | Pouch weight (g) |
| BL250 015 KIT 015 CVVH/CVVHD | Polyester/ Polypropylene Medical grade paper 60g/m ² | Polypropylene | 30.0 |
| BL250 025 KIT 025 CVVH/CVVHD | Polyester/ Polypropylene Medical grade paper 60g/m ² | Polypropylene | 30.0 |
| 10 ml Luer Lok™ BD Plastipak Syringe | Polyamied/Polyethylene Medical grade paper | N.A. | 1.3 |

| Model | Secondary packaging - Box | | | | |
|-----------------------------|--|------------------------------|--|------------------------------|---------------------------|
| | Single box | Weight ¹⁴ (kg) | Multiple box | Weight ¹⁵ (kg) | Pcs/ Box ¹⁶ |
| BL250 KIT 015 CVVH/CVVHD | White colored Rippled Cardboard 4 mm - KBM/222/B Dimensions: 560 × 377 × 71 mm | 0.9 | Avana colored Rippled Cardboard 4 mm Dimensions: 590 × 390 × 324 mm | 4.4 | 4 |
| BL250 KIT 025 CVVH/CVVHD | White colored Rippled Cardboard 4 mm - KBM/222/B Dimensions: 560 × 377 × 71 mm | 0.9 | Avana colored Rippled Cardboard 4 mm Dimensions: 590 × 390 × 324 mm | 4.5 | 4 |

¹⁴ Single box weight

¹⁵ Multiple box weight

¹⁶ Number of single boxes insert into the multiple box

Storage and disposal conditions

Storage conditions: store at temperatures between +5 and +30 degrees Celsius.

Disposal: dispose the device after treatment in accordance with applicable government and health center regulations.

Biocompatibility

Biocompatibility tests of the preassembled devices and syringe for CARPEDIEM™ machine have been performed according to EN ISO 10993-1¹⁷ and related applicable standard series.

Biocompatibility tests are the responsibility of the manufacturer of the medical devices inside the procedure packs.

¹⁷ The different ISO standard mentioned in the document are in current revision.



TECHNICAL DATA SHEET

BD Plastipak™ syringes without needles and
BD General Syringes without needle
Sterile, Single Use, Latex free

1. General Information

1.1 General

BD Plastipak™ syringe and general syringes are used for general purpose injection and aspiration of fluids from vials, ampoules and parts of the body below the surface of the skin.

Perfusion syringes, 50ml syringes, are designed for short term use in syringe pumps (active IIa devices) for the administration of pharmaceuticals. The 50 ml Catheter Tip Syringes have a long tapered tip designed to aid in irrigation or for connection to non-ISO compatible Luer connections such as nasogastric tubes.



DEAD SPACE (maximum, without needle) (except for catheter tip syringes)

| SYRINGE SIZE | 1 ml | 2ml | 5ml | 10ml | 20ml | 30ml | 50ml | 100ml |
|--------------|---------|--------|---------|--------|--------|--------|--------|--------|
| Dead Space | 0.07 ml | 0.07ml | 0.075ml | 0.10ml | 0.15ml | 0.17ml | 0.20ml | 0.20ml |



Plastipak syringes without needles TDS version Feb 2017

LUER SLIP SYRINGES

| Reference | Capacity | Description | Scale Graduation | Box (units) | Case (units) |
|-----------|----------|--------------------------------|---------------------|-------------|--------------|
| 300026 | 1 ml | Insulin 40 I.U. | International units | 100 | 800 |
| 301355 | 1 ml | Insulin 100 I.U. | International units | 100 | 800 |
| 303174 | 1 ml | Insulin 100 I.U. | International units | 120 | 960 |
| 303173 | 1 ml | Insulin 40 I.U. | International units | 120 | 960 |
| 300013 | 1 ml | Central cone | 0.01 ml | 100 | 800 |
| 303172 | 1 ml | Central cone | 0.01 ml | 120 | 960 |
| 300185 | 2/2.5 ml | Central cone | 0.1 ml | 100 | 800 |
| 302187 | 5 ml | Central cone | 0.2 ml | 100 | 400 |
| 302188 | 10 ml | Eccentric cone | 0.5 ml | 100 | 400 |
| 301183 | 20 ml | Eccentric cone | 1 ml | 60 | 240 |
| 300613 | 20 ml | Eccentric cone | 1 ml | 120 | 480 |
| 301231 | 30 ml | Eccentric cone | 1 ml | 60 | 240 |
| 300866 | 50/60 ml | Eccentric cone | 1 ml | 60 | 240 |
| 300867 | 50/60 ml | Catheter tip | 1 ml | 60 | 240 |
| 300605 | 100 ml | Catheter tip with Luer adaptor | 2 ml | 25 | 50 |
| 309654 | 60ml | Slip tip | 1 ml | 40 | 160 |

LUER LOK™ SYRINGES

| Reference | Capacity | Description | Scale Graduation | Box (units) | Case (units) |
|-----------|----------|---------------------------|------------------|-------------|--------------|
| 301189 | 20 ml | Luer Lok™ | 1 ml | 60 | 240 |
| 300629 | 20 ml | Luer Lok™ | 1 ml | 120 | 480 |
| 302830 | 20 ml | Luer Lok™ | 1 ml | 48 | 192 |
| 301229 | 30 ml | Luer Lok™ | 1 ml | 60 | 240 |
| 300865 | 50/60 ml | Luer Lok™ | 1 ml | 60 | 240 |
| 300137 | 50 ml | Luer Lok™ Perfusion | 1 ml | 50 | 100 |
| 300139 | 50 ml | Luer Lok™ Perfusion Amber | 1 ml | 50 | 100 |
| 309653 | 60 ml | Luer Lok™ | 1 ml | 40 | 160 |
| 309628 | 1 ml | Luer Lok™ | 0.01 ml | 100 | 800 |
| 309658 | 3 ml | Luer Lok™ | 0.1ml | 200 | 800 |
| 309649 | 5 ml | Luer Lok™ | 0.2ml | 125 | 500 |
| 300912 | 10 ml | Luer Lok™ | 0.2ml | 100 | 400 |
| 305959* | 10 ml | Luer Lok™ | 0.2ml | 100 | 400 |
| 300869 | 50/60 ml | Luer Lok™ Amber | 1 ml | 60 | 240 |

*305959 will be preferred to supply to European customers as this catalogue number of 10ml Luer Lok™ Plastipak is manufactured in Europe.



Plastipak syringes without needles TDS version Feb 2017

1.2 Certification

| BD REFERENCE | BD MANUFACTURER | ISO CERTIFICATION | CE MARKING | BD MANUFACTURING SITE |
|--|---|--|--|---|
| 301189, 301183, 300629, 301229, 300865, 300869, 300867, 300605, 300613, 301231, 300866, 300137, 300139 | Becton Dickinson & Company Limited Donore Road Drogheda Co. Louth Ireland | NSAI - Certificate MD 19.1609 I.S. EN ISO 13485:2012 | NSAI NB no 0050: Certificate N* 252.156 | Becton Dickinson S.A. - Camino de Valdeoliva, s/n. 28750, San Agustín del Guadalix (Madrid) Spain |
| 300026, 301355, 300013, 300185, 302187, 302188, 303172, 303173, 303174, 305959 | Becton Dickinson S.A. - Camino de Valdeoliva, s/n. 28750, San Agustín del Guadalix (Madrid) Spain | AENOR -N. ER-0264/1994 – ISO 9001:2008; AEMPS N. 2012 07 0013 EN - EN – ISO 13485:2013 | AEMPS 0318: Certificate N* 2000 06 0273 CP | Becton Dickinson S.A. - Camino de Valdeoliva, s/n. 28750, San Agustín del Guadalix (Madrid) Spain |
| 309628*, 309658, 309649, 300910, 300911, 300912, | Becton, Dickinson and Company 1 Becton Drive Franklin Lakes, NJ 07417, USA | NSAI - ISO 9001 :20008 Certificate MD19.2305 NSAI ISO 13485 :2012 Certificate MD19.2305 | NSAI 0050: Certificate N* 252.231 | Becton, Dickinson and Company Route 7 & Grace Way, Canaan CT 06018 USA |
| 309653, 309654, 302830 | Becton, Dickinson and Company 1 Becton Drive Franklin Lakes, NJ 07417, USA | NSAI - ISO 9001 :20008 Certificate MD19.2305 NSAI ISO 13485 :2012 Certificate MD19.2305 | NSAI 0050: Certificate N* 252.231 | BD Medical - Medical Surgical Systems 2153 12th Avenue Columbus, NE 68602 USA |

*Catalogue number 309628 used to be manufactured in BD Singapore Branch, 30 Tuas Avenue 2, Singapore 639461. No changes to form, fit or function when transferred to BD Canada.



Plastipak syringes without needles TDS version Feb 2017

1.3 Material

| COMPONENT | MATERIAL |
|-----------------------|--|
| SYRINGE | |
| Barrels, plunger rods | POLYPROPYLENE |
| Barrel cat# 309628 | POLYCARBONATE |
| Stoppers | LATEX FREE ELASTOMER |
| Lubricant | MEDICAL GRADE SILICONE OIL, <0.25 mg/cm ² |

BD Plastipak™ amber syringes, such as 300139 and 300869, have the barrel colored to reduce U.V. light for administration of light sensitive medications. The light transmission has been characterized as per the transparency test (method 1) as described in the Japanese Pharmacopeia XVI.

According to such method the light transmissibility (%) is characterized under a UV light source emitting at 450nm. The light transmission is 5.6± 0.2% (mean ± standard deviation) according to the transparency test.

1.4 Material of concern

Materials of concern are chemicals or substances that have been identified as having the potential to cause long term effects on humans or the environment.

| MATERIAL | COMMENT |
|-------------------------------------|---|
| Phthalates | The products do not contain phthalates. No DEHP, CAS number 117-81-1, EC number 204-211-0, intentionally added |
| Latex | The products do not contain natural latex. |
| Bisphenol A | The products do not contain Bisphenol A. Catalogue number 309628 contain polycarbonate and hence Bisphenol A |
| Substances of animal origin BSE/TSE | The products were assessed for TSE (Transmissible Spongiform Encephalopathy) contamination risk. The raw materials used in the manufacture of this device do not contain any animal tissue but may contain very small amounts of animal derived raw materials. This product is manufactured using polymer resins which may contain very small amounts of surfactants or fatty acids derived from tallow. Our resin suppliers have confirmed that these tallow derived materials have been produced using multiple cycles of conditions at least as rigorous (and normally more rigorous) as those specified in Annex C.5 of EN ISO 22442-1. Therefore, the raw materials meet or exceed the requirements of EN ISO 22442-1. Therefore, the raw materials meet or do not to present any risk with respect to TSE or other animal borne disease |
| Polyvinyl chloride (PVC) | The products do not contain polyvinyl chloride |

1.5 REACH information

BD maintains an active REACH compliance program and works closely with its supply base on an ongoing basis with a view to obtaining information on REACH Substances of Very High Concern (“SVHC”) through regular communication and exchange



1.6 Biocompatibility

BD Medical products comply with the requirements of the standard for biocompatibility of medical devices, ISO 10993-1 Biological Evaluation of Medical Devices - Part 1: Evaluation and Testing.

1.7 Sterilization

- **Ethylene Oxide Sterilization** following EN ISO 11135-1. ETO residues are within applicable regulations. All references except references below are sterilized with EO
- **Radiation Sterilization** following EN ISO 11137-1 References sterilized with radiation: 309628, 309658, 309649, 300910, 300911, 300912, 302830, 309653 and 309654.

1.8 Shelf life

Shelf life 5 years for all catalogue numbers except 300605.
Catalogue number 300605, 100ml Catheter tip, has a shelf life of 18 months.
No special storage or transportation condition. Recommendations are to store in room temperature, in dry and warm place and not exposed to strong light.

1.9 Standards

| HARMONISED STANDARDS | |
|-----------------------------|--|
| EN 556-1:2001/ AC:2006 | Sterilisation of Medical Devices – requirements for medical devices to be labelled "sterile". |
| EN 980: 2008 | Graphical Symbols for use in the labelling of medical devices. |
| BS EN 1041+A1: 2013 | Terminology, symbols and information provided with medical devices. Information supplied by the manufacturer with medical devices |
| EN 1707:1996 | Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment - Lock fittings |
| EN 20594-1:1993/AC:1996 | Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment - Part 1: General requirements |
| EN ISO10993-series | Biological evaluation of medical devices |
| EN ISO 11135-1:2007 | Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices |
| EN ISO 11137-1 | Sterilization of health care products - Radiation. Part1. Requirements for development, validation and routine control of sterilization process for medical devices |
| EN ISO 11137-2 | Sterilization of health care products – Radiation. Part2. Establishing the sterilization dose |
| EN ISO 11138-2:2009 | Sterilization of health care products - Biological Indicators - Part 2: Biological indicators for ethylene oxide sterilization processes |
| EN ISO 11607-1:2009 | Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems |
| EN ISO 11607-2:2006 | Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes |
| EN ISO 11737-1:2006/AC:2009 | Sterilization of medical devices - Microbiological methods - Part 1: Determination of a population of microorganisms on products |

This document is approved electronically
This document can be changed without further notification



Plastipak syringes without needles TDS version Feb 2017

| HARMONISED STANDARDS, continue | |
|--------------------------------|---|
| EN ISO 11737-2:2009 | Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process (ISO 11737-2:2009) |
| EN ISO 13485:2012/AC:2012 | Medical devices – Quality management Systems Requirements for Regulatory Purposes |
| EN ISO 14155:2011 | Clinical investigation of medical devices for human subjects - Good clinical practice |
| EN ISO 14971:2012 | Medical Devices. Application of risk management to medical devices |

| NON HARMONISED STANDARD | |
|-------------------------|---|
| IS EN ISO 7864-1: 1996 | Sterile hypodermic needles For Single Use |
| IS EN ISO 7886-1:1998 | Sterile hypodermic syringes For Single use - Part 1:Syringes for manual use <i>See Note 1 Below</i> |
| EN ISO 7886-2:1998 | Sterile Hypodermic Syringes for Single Use. Part 2: Syringes for Use with Power-Driven Syringe Pumps. <i>See notes 2 and 3 below</i> |
| ISO 594-1:1993 | Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment - Part 1: General requirements |
| ISO 594-2:1998 | Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment – Part 2: Lock fittings |
| ISO 9626: 1995 | Stainless steel needle tubing for the manufacture of medical devices |
| ISO 13485:2003 | Medical devices – Quality management Systems Requirements for Regulatory Purposes |
| ISO 14644-1:1999 | Cleanrooms and associated controlled environments Part 1: Classification of air cleanliness |
| ISO 15223-1:2012 | Sterilisation of Medical Devices – requirements for medical devices to be labelled “sterile”. |
| ISO 10993-2:2009 | Biological Evaluation of Medical Devices Part 2 |
| ISO10993-10:2009 | Biological evaluation of medical devices - Part 10: Tests for Irritation and delayed-type hypersensitivity |
| ISO 2859-1:1999 | Sampling procedures for inspection by attributes – Part 1: Sampling schemes Indexed by acceptance quality limit (AQL) for lot-by-lot inspection |

Notes:

- Subclause 14.1, 50ml Perfusion is non-compliant; Dead space approx. double the standard requirement (0.38ml versus 0.2ml maximum)
- Subclause 14, Plastipak syringes (including Perfusion) do not comply as in our opinion the requirements in ISO 7886-2 section 14 are most likely not applicable to the current pump landscape and clinical requirements in the market place. Also, some of the requirements in section 14 appear questionable from a technical perspective, based on the latest state of BD’s knowledge about syringe performance on pumps and related testing capabilities
- Subclauses 16.1, following statement not on the labelling for Plastipak Luer Lok “ Suitable for use with power-driven syringe pumps or equivalent
INSULIN GRADUATED SYRINGE ALSO MEETS ISO 8537 Sterile single-use syringes, with or without needle, for insulin



Plastipak syringes without needles TDS version Feb 2017

1.10 Classification

- **Class I:** 300026, 301355, 300013, 300185, 302187, 302188, 303172, 303173, 303174 and 305959, Rule 2, Annex IX of Medical Devices Directive 93/42/EEC as amended.
- **Class I:** 309628, 309658, 309649, 300910, 300911, 300912, 309653, 309654 and 302830, Rule 1, Annex IX, Section III of the Medical Device Directive 93/42/EEC as amended
- **Class IIa:** 301189, 301183, 300629, 301229 300865, 300869, 300867, 300605, 300613, 301231, 300866, 300137 and 300139, Rule 2, Annex V and VII of the Medical Devices Directive 93/42/EEC as amended.

1.11 GMDN code

GMDN code 47017: General purpose syringes.

1.12 Good Manufacturing practices

The entire manufacturing and testing processes are following the Good Manufacturing Practices as specified below

- Incoming raw materials are verified via material inspection and testing and our suppliers are approved via our vendor management system.
- In addition to the automatic on-line inspections, in-process inspections are performed in addition to final product testing to ensure compliance with approved specifications.
- The manufacturing and testing details of each batch of product are recorded on a batch record which is retained in accordance with our document control procedures
- BD operates a system of Internal and external audits to maintain compliance
- BD confirms that it will continue to adhere to relevant international standards in designing and manufacturing its products.
- BD reserves the right to use the internal change control procedure to change raw material suppliers and production process

1.13 Others

- The EU representative, for syringes which BD Manufacturer is BD Franklin Lakes such as 309628, 309658, 309649, 300910, 302830, 300911, 300912, 309653 and 309654 is Becton Dickinson Distribution Center, Laagstraat 57, B-9140 Temse -Belgium. Other syringes are produced by a European manufacturer.
- (Material) Safety Data Sheets are not required for this product
- Certificate of Food Contact (*COMMISSION REGULATION (EU) No. 10/2011 of January 14th, 2011 concerning materials and plastic objects intended to get in touch with foodstuffs*) is not required as BD products are used for general purpose injection and aspiration of fluids from vials, ampoules and parts of the body below the surface of the skin.

Becton Dickinson and Company
BD Medical,
Medication and Procedural Solution

BD Switzerland Sàrl
Terre Bonne Park - A4
Route de Crassier 17
1262 Eysins,
Switzerland



bd.com

Plastipak syringes without needles TDS version Feb 2017

2. Packaging

2.1 Packaging material

| | |
|---------------|---|
| PACKAGING | |
| Web packaging | Polyamied/polyethylene, Medical grade paper |
| Ink | Printing Ink |
| Box | Hard Paper |

LABELS: according to European Medical Device directive, multilingual

2.2 Example labeling

Legal Manufacturer and manufacturing site: San Agustin del Guadalix
Example Unit pack cat.no 301355, from document DGW1318

Legal Manufacturer: Drogheda and manufacturing site San Agustin del Guadalix
Example Unit pack cat.no 300867, from document DGW1086

ment is approved electronically
ment can be changed without further notification

Becton Dickinson and Company
BD Medical,
Medication and Procedural Solution

BD Switzerland Sàrl
Terre Bonne Park - A4
Route de Crassier 17
1262 Eysins,
Switzerland

bd.com



Plastipak syringes without needles TDS version Feb 2017

Legal Manufacturer: Franklin Lakes and Manufacturing site Cnaan
Example Unit pack cat.no 309628, from document number DGW757

1ml Luer-Lok™ Syringe

BD 1ml Syringe
Luer-Lok™ Tip

Jeringa • Seringa • Seringue • Spritze •
Siringa • Spuit • Spruta • Sprøjte •
Ruisku • Σύριγγα • Sprøyte •
Strzykawka • Injekcijska brizgalka •
Injekčná striekačka • Süstal •
Fecskendő • Švirkštas • Stříkačka •
Şırıce • Şiringa • Шприц • Štrcaljka •
Seringă • Спринцовка • Шприц



REF 309628



0050

STERILE R STERILIZED USING
IRRADIATION

 DO NOT
REUSE

 CAUTION, CONSULT
ACCOMPANYING DOCUMENTS



(01)00382903096282

MANUFACTURER
Becton, Dickinson and Company, 1 Becton Drive, Franklin Lakes, NJ 07417 USA
Made in USA, www.bd.com BD, BD Logo and BD Luer-Lok are
trademarks of Becton, Dickinson and Company, © 2016 BD

DGW75704
500011886

This document is approved electronically
This document can be changed without further notification

Page 9 of 9

Advancing the world of health

Medtronic

Bellco S.r.l.
Via Camurana, 1
41037 Mirandola (MO) Italy
Tel: +39 0535 29111

www.bellco.net

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.
© 2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.
™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 22-bl250-kit-015-and-025-cvvh-cvvh-d-data-sheet-en-we-7727216