

Mahurkar™ Chronic Carbothane™ Catheter


Reliability and performance.

Performance

- 14.5Fr Double-D™ lumen design maximizes the cross-sectional area within the lumen
- Silicone extensions allow visibility of blood
- Ultem™ adapters withstand repeated use

Safety

- VenaTrac™ over-the-wire insertion stylets reduce risk of air embolism and blood loss¹
- Safety tunneler allows for ease of insertion
- Large laser-printed volumes for easy reading



VenaTrac™ stylets:

- Allow for sheathless catheter insertions¹
- Promote over-the-wire catheter exchanges¹



Valved pull-apart safety sheath

Ordering information

		Priming Volumes			
	Code	Description	Arterial (Red)	Venous (Blue)	Qty Case
Mahurkar™ Chronic Carbothane™ Catheter					
Sport Pack	8888119385	19/36 with sideholes	1.6 mL	1.7 mL	5
	8888123435	23/40 with sideholes	1.8 mL	1.9 mL	5
	8888128485	28/45 with sideholes	2.0 mL	2.1 mL	5
	8888133535	33/50 with sideholes	2.2 mL	2.3 mL	5
Kits	8888145251	19/36 with sideholes	1.6 mL	1.7 mL	5
	8888145252	23/40 with sideholes	1.8 mL	1.9 mL	5
	8888145253	28/45 with sideholes	2.0 mL	2.1 mL	5
	8888145254	33/50 with sideholes	2.2 mL	2.3 mL	5
Mahurkar™ Chronic Carbothane™ Catheter with VenaTrac™ Over-the-Wire Insertion Stylets (One Pair)					
Sport Pack	8888119362	19/36 with sideholes	1.6 mL	1.7 mL	5
	8888123402	23/40 with sideholes	1.8 mL	1.9 mL	5
	8888128452	28/45 with sideholes	2.0 mL	2.1 mL	5
	8888133502	33/50 with sideholes	2.2 mL	2.3 mL	5
Kits	8888145027	19/36 with sideholes	1.6 mL	1.7 mL	5
	8888145028	23/40 with sideholes	1.8 mL	1.9 mL	5
	8888145029	28/45 with sideholes	2.0 mL	2.1 mL	5
	8888145030	33/50 with sideholes	2.2 mL	2.3 mL	5

Components	Qty	Sport pack	Kit
Mahurkar™ Chronic Carbothane™ Catheter	1	•	•
16Fr valved pull-apart safety sheath	1	•	•
Tunneler, 21 cm	1	•	•
12Fr tissue dilator	1	•	•
14Fr tissue dilator	1	•	•
Sealing caps	2	•	•
18 G x 7 cm introducer needle	1		•
J/Straight .038"x 120 cm guidewire	1		•
No. 11 scalpel	1		•
Telfa™ Island Dressings	2		•
4" x 4" cotton gauze sponges	4		•
12 mL syringe	1		•

For more information, visit [MozarcMedical.com](https://www.MozarcMedical.com)

1. Mojibian HR, Oneal P, Eliseo D, Tal MG. VenaTrac™ device for over-the-wire placement of chronic dialysis catheters. *J Vasc Access.* 2010;11(4):342-345.

Mahurkar™ Chronic Carbothane Catheter

The Mahurkar 14.5 Fr Chronic Carbothane catheter is intended for acute and chronic hemodialysis, apheresis, and infusion. It may be inserted either percutaneously or by cutdown.

Do not use the catheter in thrombosed vessels or for subclavian puncture when ventilator is in use.

See the device manual for detailed information regarding the implant procedure, indications, contraindications, warnings, precautions, and potential complications/adverse events. For further information, please call Mozarc Medical, in the United States call 1-800-962-9888 and for international inquiries, call 508-261-8000, and/or consult www.MozarcMedical.com.

Copyright © 2023 Mozarc Medical Holding LLC. Mozarc, Mozarc Medical, the Mozarc Medical logos, and Empowering Patients. Enriching Lives. are trademarks of Mozarc Medical. MAHURKAR™, Double-D™, and the Double-D design are U.S. registered trademarks of Sakharam D. Mahurkar, used under license. ™ Third-party brands are trademarks of their respective owners. 08/2023 – US-RC-2300052 – [WF#10213676]

710 Medtronic Parkway
Minneapolis, MN 55432
800.962.9888
508.261.8000

[MozarcMedical.com](https://www.MozarcMedical.com)

Mozarc
Empowering patients.
Enriching lives. **medical**