


WHY MANAGE THE RISK WHEN WE CAN REDUCE IT?




16 Fr valved pull-apart introducer sheath

Palindrome™ Precision dialysis catheter kits



DECREASE THE RISK OF AIR EMBOLISM

The widespread use of valved peel-away sheaths prevents air embolism when the catheter is placed.¹ The valve contained inside the 16 Fr pull-apart introducer sheath automatically activates after removing the dilator, without the need of any additional maneuver from the operator.



INTUITIVE SECUREMENT


The locking mechanism secures in place the dilator to the 16 Fr pull-apart introducer sheath for a safe introduction by simply rotating the dilator.



EASY TO USE

The ergonomic tabs include ribs to provide grip and easily pull apart the sheath while inserting the dialysis catheter.

The tapered tip of the 16 Fr pull-apart introducer sheath allows for a smooth transition to the dilator and eases the introduction of the assembled device into the skin and vessel of the patient.




Medtronic

Palindrome™ Precision dialysis catheter kits including the new 16 Fr pull apart introducer sheath

Product codes	Description
8888145014P	14.5 Fr Palindrome™ Precision catheter Kit, 19 cm tip to cuff, 36 cm
8888145015P	14.5 Fr Palindrome™ Precision catheter Kit, 23 cm tip to cuff, 40 cm
8888145016P	14.5 Fr Palindrome™ Precision catheter Kit, 28 cm tip to cuff, 45 cm
8888145017P	14.5 Fr Palindrome™ Precision catheter Kit, 33 cm tip to cuff, 50 cm
8888146044P	14.5 Fr Palindrome™ Precision catheter Kit, 44 cm tip to cuff, 61 cm
8888145018P	14.5 Fr Palindrome™ Precision catheter Kit, 55 cm tip to cuff, 72 cm
8888145039P	14.5 Fr Palindrome™ Precision catheter Kit with VenaTrac, 19 cm tip to cuff, 36 cm
8888145040P	14.5 Fr Palindrome™ Precision catheter Kit with VenaTrac, 23 cm tip to cuff, 40 cm
8888145041P	14.5 Fr Palindrome™ Precision catheter Kit with VenaTrac, 28 cm tip to cuff, 45 cm
8888145042P	14.5 Fr Palindrome™ Precision catheter Kit with VenaTrac, 33 cm tip to cuff, 50 cm
8888145043CP	14.5 Fr Palindrome™ Precision Heparin-Coated catheter Kit, 19 cm tip to cuff, 36 cm
8888145044CP	14.5 Fr Palindrome™ Precision Heparin-Coated catheter Kit, 23 cm tip to cuff, 40 cm
8888145045CP	14.5 Fr Palindrome™ Precision Heparin-Coated catheter Kit, 28 cm tip to cuff, 45 cm
8888145046CP	14.5 Fr Palindrome™ Precision Heparin-Coated catheter Kit, 33 cm tip to cuff, 50 cm
8888145057CP*	14.5 Fr Palindrome™ Precision Heparin-Coated and Silver ion sleeve catheter Kit, 19 cm tip to cuff, 36 cm
8888145048CP*	14.5 Fr Palindrome™ Precision Heparin-Coated and Silver ion sleeve catheter Kit, 23 cm tip to cuff, 40 cm
8888145049CP*	14.5 Fr Palindrome™ Precision Heparin-Coated and Silver ion sleeve catheter Kit, 28 cm tip to cuff, 45 cm
8888145050CP*	14.5 Fr Palindrome™ Precision Heparin-Coated and Silver ion sleeve catheter Kit, 33 cm tip to cuff, 50 cm
8888145043P*	14.5 Fr Palindrome™ Precision Heparin-Coated catheter Kit, 19 cm tip to cuff, 36 cm
8888145044P*	14.5 Fr Palindrome™ Precision Heparin-Coated catheter Kit, 23 cm tip to cuff, 40 cm
8888145045P*	14.5 Fr Palindrome™ Precision Heparin-Coated catheter Kit, 28 cm tip to cuff, 45 cm
8888145046P*	14.5 Fr Palindrome™ Precision Heparin-Coated catheter Kit, 33 cm tip to cuff, 50 cm
8888145057P*	14.5 Fr Palindrome™ Precision Heparin-Coated and Silver ion sleeve catheter Kit, 19 cm tip to cuff, 36 cm
8888145048P*	14.5 Fr Palindrome™ Precision Heparin-Coated and Silver ion sleeve catheter Kit, 23 cm tip to cuff, 40 cm
8888145049P*	14.5 Fr Palindrome™ Precision Heparin-Coated and Silver ion sleeve catheter Kit, 28 cm tip to cuff, 45 cm
8888145050P*	14.5 Fr Palindrome™ Precision Heparin-Coated and Silver ion sleeve catheter Kit, 33 cm tip to cuff, 50 cm
8888145062P*	14.5 Fr Palindrome™ Precision with Silver ion sleeve catheter Kit, 19 cm tip to cuff, 36 cm
8888145063P*	14.5 Fr Palindrome™ Precision with Silver ion sleeve catheter Kit, 23 cm tip to cuff, 40 cm
8888145064P*	14.5 Fr Palindrome™ Precision with Silver ion sleeve catheter Kit, 28 cm tip to cuff, 45 cm
8888145065P*	14.5 Fr Palindrome™ Precision with Silver ion sleeve catheter Kit, 33 cm tip to cuff, 50 cm
8888146144P*	14.5 Fr Palindrome™ Precision with Silver ion sleeve catheter Kit, 44 cm tip to cuff, 61 cm
8888145066P*	14.5 Fr Palindrome™ Precision with Silver ion sleeve catheter Kit, 55 cm tip to cuff, 72 cm

* Product codes not available in Europe.

1. Peter R. Bream Jr.: Update on Insertion and Complications of Central Venous Catheters for Hemodialysis. <http://dx.doi.org/10.1055/s-0036-1572547>. ISSN 0739-9529.



Medtronic

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2023 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. EMEA-RC-2100003 - 23-palindrome-precision-16-fr-pull-apart-brochure - 4893483