

# Bridge the gap to peritoneal dialysis care options

## Ritus™ peritoneal dialysis tunneling tool kit

### Expand access

Enable the external portion of a peritoneal dialysis (PD) catheter to be temporarily buried within subcutaneous abdominal tissue — in advance of needing peritoneal dialysis therapy.

For physicians looking to expand the adoption of home dialysis therapy, specifically PD, the Ritus™ PD tunneling tool kit provides the components needed for proactive catheter insertion.


The Ritus™ tunneling tool kit includes:

- Stainless steel tunneling device
- Removable plastic cap (for catheter attachment to pull through the tunnel track)
- Titanium catheter plug (to seal the catheter end)

### Remove barriers

Embedding PD catheters in patients with chronic kidney disease (CKD) can allow for:

- More efficient surgical scheduling<sup>1</sup>
- Greater patient acceptance<sup>1</sup>
- Avoid the need for the initiation of hemodialysis with a catheter.<sup>1</sup>



# Ordering information

For full portfolio details, please contact your local Mozarc Medical sales representative.

Visit [medtronic.com/PeritonealDialysis](https://www.medtronic.com/PeritonealDialysis)



Item number	Description	Quantity
8817791002	Ritus™ tunneling tool kit	1/pack

For patients who are future peritoneal dialysis treatment candidates, the Ritus™ tunneling tool kit allows the operator to temporarily bury the external portion of the PD catheter in advance of dialysis.

Explore options for compatible peritoneal dialysis catheters with the Argyle™\* peritoneal dialysis family.

## Related PD catheters and kits:

Item number	Description	Quantity
8811313015	Peritoneal Curl catheter, two cuffs, 57 cm	1/case
8817278007	Peritoneal Curl catheter kit, two cuffs, 57 cm	5/case
8811313010	Peritoneal Curl catheter, two cuffs, 62 cm	1/case
8817278006	Peritoneal Curl catheter kit, two cuffs, 62 cm	5/case
8888413823	Peritoneal Swan Neck Curl catheter kit, two cuffs, left 62.5 cm	1/case
8888413831	Peritoneal Swan Neck Curl catheter kit, two cuffs, right 62.5 cm	1/case
8888111132	Peritoneal Swan Neck pre-sternal two cuffs, Curl Cath one cuff tray	1/case

Scan to view a full list of PD catheter configuration options



## Related PD insertion accessories:

Item number	Description	Quantity
8888412305	Ritus™ preoperative peritoneal dialysis catheter stencil	5/pack
8888415679	Faller tunneling stylet	1/case
8888415661	Insertion straightening stylet 61 cm	1/case
8888415604	Titanium two-part adapter	1/case

### Ritus™ peritoneal dialysis catheter tunneling tool kit

The Ritus™ tunneling tool kit is used to both create a tunnel tract through subcutaneous tissue and seal the end of the catheter. This allows the external portion of the Covidien™ Argyle™ peritoneal dialysis catheter to be temporarily buried in advance of dialysis for patients that are candidates for delayed peritoneal dialysis treatment. Risks may include peritonitis and exit site infection.

### Contraindications

Do not embed the external portion of the catheter in a separate procedure than the catheter implantation. Do not use this product if the patient is not a candidate for delayed peritoneal dialysis treatment or requires immediate peritoneal dialysis treatment. Not for use in pediatric, infant and neonate.

1. Crabtree JH, Burchette RJ, Siddiqi RA. Embedded Catheters: Minimizing Excessive Embedment Time and Futile Placement while Maintaining Procedure Benefits. *Perit Dial Int.* 2015;35(5):545-551. doi:10.3747/pdi.2013.00301

© 2023 Mozarc Medical Holding LLC. Mozarc, Mozarc Medical, the Mozarc Medical logos, and Empowering Patients. Enriching Lives. are trademarks of Mozarc Medical. Covidien and Medtronic are trademarks of Medtronic used under license.™\* Third-party brands are trademarks of their respective owners. 08/2023 – US-RC-2300046 – [WF#10010561]

Medtronic Parkway 710  
Minneapolis, MN 55432  
[MozarcMedical.com](https://www.MozarcMedical.com)

**Mozarc**  
Empowering patients. Enriching lives. **medical**