


Argyle™ peritoneal dialysis catheters (PD)

# Explore PD with a personalized fit

## Tenckhoff feature



- Straight extension
- (1-2) Felt cuffs



12 different size/configurations

## Curl cath feature



- Curled mass
- Inflow-outflow diffusion holes
- (1-2) Felt cuffs



13 different size/configurations

## Swan neck feature



- Pre-formed angle
- Inflow-outflow diffusion holes
- (1-2) Felt cuffs
- Presternal option



14 different size/configurations



## Missouri feature

- (1) Silicone bead
- (1) Felt disc (internal; rectus muscle)
- (2) Felt cuffs with secondary subcutaneous felt cuff (exit site seal)





## Oreopoulos-Zellerman feature

- (2) Two silicone discs (internal)



## Moncrief- Popovich feature

- Enlarged subcutaneous felt cuff
- Supports embedded insertion technique



**Learn more:** See the full brochure for all size, length, and cuff options

Tenckhoff Configurations:	Length (cm)
One attached cuff (subcutaneous or preperitoneal position)	41, 42 , 46
Two attached cuffs (one subcutaneous and preperitoneal position)	42, 47
Two disks with one or two cuffs - Oreopolulos-Zellerman Missouri	41

Swan neck Style/ Configuration	Length (cm)	Placement
Swan Neck Tenckhoff, 2 cuff	43	left
Swan Neck Curl Cath, 2 cuff	62.5	Right and left side
Swan Neck Oreopoulos-Zellerman Catheters, 2 cuff	44.5	Right and left side
Swan Neck Missouri Catheters, 2 cuff	44.5	Right and left side
Swan Neck Curl Cath Missouri Catheter, 2 Cuff	62.5	Right and left side
Moncrief-Popovich Swan Neck Curl Cath Catheters, 2 cuff	62.5	
Swan Neck Presternal Catheters 2 cuff	60.3 (upper segment) 52 Cm (Lower segment); 112.3 cm total	

Curl Cath peritoneal dialysis catheter	
Style/Configuration	Lengths (cm)
1 Cuff	57, 60, 62
2 Cuff	57, 62

**Lumen diameter:** 15Fr/Ch  
**Material:** Translucent silicone rubber

Pediatric peritoneal dialysis catheter	
Style	Lengths (cm)
Tenckhoff	30, 31, 37
Tenckhoff Curl Cath	39
Swan Neck Tenckhoff	37.5
Swan Neck Curl Cath	38.9, 42, 43, 59,
Swan Neck Missouri	38
Swan Neck Curl Cath Missouri	42

©2024 Mozarc Medical Holding LLC. Mozarc, Mozarc Medical, the Mozarc Medical logos, and Empowering Patients. Enriching Lives. are trademarks of Mozarc Medical. Covidien and Medtronic are trademarks of Medtronic used under license.

™\* Third-party brands are trademarks of their respective owners. 02/2024 – US-RC-2400039 [IH]

800.962.9888  
508.261.8000

710 Medtronic Parkway  
Minneapolis, MN 55432

See PD catheter and accessories details at:  
[Medtronic.com/dialysis-access](https://www.Medtronic.com/dialysis-access)  
[MozarcMedical.com](https://www.MozarcMedical.com)

