


# The range needed for the right peritoneal access

Explore the Argyle™\*  
Peritoneal Dialysis  
Catheter portfolio



Different patients have different anatomies and needs, so why settle for a less than ideal fit? With Argyle™\* peritoneal dialysis catheters, you don't have to. Our variety of configurations and lengths gives you options designed to support long-term access preservation and patency.

| Available Peritoneal Components |                                                                                               |                                                                                                         |
|---------------------------------|-----------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|
| Internal end                    | <br>Straight | <br>Coiled           |
| Fixation components             | <br>1 cuff   | <br>2 cuff           |
| Transmural segment              | <br>Straight | <br>Pre-formed bend |



1. Radiopaque Silicone  
Identify under x-ray.

2. Felt Cuff

3. Bend  
Swan neck options.

4. Missouri Bead + Felt Cuff

5. Oreopoulos-Zellerman  
Option with two  
silicon disks.

6. Straight (Tenckhoff)  
Indwelling portion  
of catheter.

7. Curled  
Indwelling portion  
of catheter.

8. Moncrief-Popovich

# Our complete peritoneal access portfolio

Peritoneal Access products can be ordered in either a catheter configuration or a catheter kit configuration.

## Catheter Configuration Includes:

Sterile contents included in the kit

- (1) Peritoneal Catheter
- (1) Beta-Cap Adapter
- (1) Beta-Cap Clamp

When ordering a single catheter, you may also purchase these accessories:

## SINGLE CATHETER ACCESSORIES:


### 16 Fr Pull-Apart Introducer Set

| Code                                                             | Includes                            | Ship Case |
|------------------------------------------------------------------|-------------------------------------|-----------|
| 8815544010<br><br>† Made with PolyTetra<br>FluoroEthylene (PTFE) | (1) 16 Fr Pull-Apart Sheath/Dilator | 5         |
|                                                                  | (1) 18 ga Introducer Needle         |           |
|                                                                  | (1) 12 mL Syringe                   |           |
|                                                                  | (1) J/Straight 0.038" Guide Wire    |           |




### Tunneling Stylet, sterile, single use (choose one tunneling stylet per catheter)

| Code       | Description                            | Ship Case |
|------------|----------------------------------------|-----------|
| 8817791001 | Atraumatic Peritoneal Tunneling Stylet | 5         |
| 8888415679 | Faller Tunneling Stylet                | 1         |



### Straightening Stylet, Non-sterile, Autoclavable

| Code       | Description                           | Ship Case |
|------------|---------------------------------------|-----------|
| 8888415661 | Insertion/Straightening Stylet, 61 cm | 1         |



## Catheter Kit Configuration Includes:

### Sterile contents included in the kit

- (1) Peritoneal Catheter
- (1) Beta-Cap Adapter, Cap and Clamp
- (1) 16 Fr Pull-Apart Sheath†Dilator
- (1) 18 ga Introducer Needle
- (1) J/Straight Guide Wire
- (1) 12 mL Syringe
- (1) #11 Scalpel
- (1) Atraumatic Peritoneal Tunneling Stylet
- (1) 4" × 4" Gauze Sponges

† Made with PolyTetra FluoroEthylene (PTFE)

When ordering a catheter kit, you may also purchase this accessory:

## Catheter Kit Accessory:

### Straightening Style, Non-sterile, Autoclavable (choose one stylet per catheter)

| Code       | Includes                              | Ship Case |
|------------|---------------------------------------|-----------|
| 8888415661 | Insertion/Straightening Stylet, 61 cm | 1         |



# Tenckhoff Peritoneal Catheters



## Argyle™\* Tenckhoff Peritoneal Catheters

Choose from a variety of catheter lengths and cuff placements and types, or choose the Argyle™\* Universal Tenckhoff catheter – a catheter that can be trimmed to fit most patient sizes. The Tenckhoff catheter is offered in the following cuff configurations:


- One attached cuff – in either the subcutaneous or preperitoneal position
- Two attached cuffs – one subcutaneous and one preperitoneal
- Two disks with one or two cuffs – Oreopoulos-Zellerman, Missouri catheter

### Tenckhoff Catheters - One Cuff


| Type     | Code       | Description                                     | Ship Case |
|----------|------------|-------------------------------------------------|-----------|
| Catheter | 8814843001 | Tenckhoff Catheter, 1 Preperitoneal Cuff, 41 cm | 1         |



| Type     | Code       | Description                                    | Ship Case |
|----------|------------|------------------------------------------------|-----------|
| Catheter | 8810889003 | Tenckhoff Catheter, 1 Subcutaneous Cuff, 42 cm | 1         |





| Type     | Code       | Description                                                    | Ship Case |
|----------|------------|----------------------------------------------------------------|-----------|
| Catheter | 8814843002 | Tenckhoff Catheter, Universal, 1 Preperitoneal Cuff, 46 cm     | 1         |
| Kit      | 8817278001 | Tenckhoff Catheter, Universal, 1 Preperitoneal Cuff Kit, 46 cm | 5         |




Tenckhoff Catheters - Two Cuffs


| Type     | Code       | Description                        | Ship Case |
|----------|------------|------------------------------------|-----------|
| Catheter | 8888411405 | Tenckhoff Catheter, 2 Cuffs, 41 cm | 1         |




| Type     | Code       | Description                            | Ship Case |
|----------|------------|----------------------------------------|-----------|
| Catheter | 8810888003 | Tenckhoff Catheter, 2 Cuffs, 42 cm     | 1         |
| Kit      | 8817278008 | Tenckhoff Catheter, 2 Cuffs Kit, 42 cm | 5         |



| Type     | Code       | Description                                   | Ship case |
|----------|------------|-----------------------------------------------|-----------|
| Catheter | 8810888012 | Tenckhoff Catheter, Universal, 2 Cuffs, 47 cm | 1         |
| Catheter | 8888411447 | Tenckhoff Catheter, 2 Cuffs, 47 cm            | 1         |





| Type     | Code       | Description                               | Ship Case |
|----------|------------|-------------------------------------------|-----------|
| Catheter | 8810887003 | Tenckhoff Catheter, 2 Mobile Cuffs, 42 cm | 1         |



Oreopoulos-Zellerman Missouri Catheters

| Type     | Code       | Description                                                                      | Ship Case |
|----------|------------|----------------------------------------------------------------------------------|-----------|
| Catheter | 8888411009 | Oreopoulos-Zellerman Catheter, Missouri Catheter 1 Bead, 2 Disks, 2 Cuffs, 41 cm | 1         |



# Swan Neck Peritoneal Catheters

## Argyle™\* Swan Neck Peritoneal Catheters



Choose from a variety of catheter options: Tenckhoff, Curl Cath, Missouri, Oreopoulos-Zellerman and Presternal. Paediatric catheters are also available.

- The Swan Neck catheter employs a pre-formed angle.
- The felt cuffs anchor the catheter subcutaneously.
- A radiopaque stripe on the front of each catheter makes it visible to differentiate under x-ray.
- Stencils are available for skin markings prior to surgery. Each stencil reflects the exact shape of the catheter.


## Swan Neck Tenckhoff Catheters

Swan Neck Tenckhoff catheter is a standard Tenckhoff with a pre-formed angle.

| Type     | Code       | Description                                        | Ship Case |
|----------|------------|----------------------------------------------------|-----------|
| Catheter | 8888412007 | Swan Neck Tenckhoff Catheter, 2 Cuffs, Left, 43 cm | 1         |
| Catheter | 8888425728 | Swan Neck Tenckhoff Catheter, 2 Cuffs, Left, 48 cm | 1         |




| Type     | Code       | Description                                         | Ship Case |
|----------|------------|-----------------------------------------------------|-----------|
| Catheter | 8888412015 | Swan Neck Tenckhoff Catheter, 2 Cuffs, Right, 43 cm | 1         |




Swan Neck Curl Cath Catheters

Swan Neck Curl Cath catheter employs the weight of the curl and the pre-formed bend. The outside portion of the catheter is positioned downward.

| Type     | Code       | Description                                              | Ship Case |
|----------|------------|----------------------------------------------------------|-----------|
| Catheter | 8888413807 | Swan Neck Curl Cath Catheter, 2 Cuffs, Left, 62.5 cm     | 1         |
| Kit      | 8888413823 | Swan Neck Curl Cath Catheter Kit, 2 Cuffs, Left, 62.5 cm | 1         |





| Type     | Code       | Description                                               | Ship Case |
|----------|------------|-----------------------------------------------------------|-----------|
| Catheter | 8888413815 | Swan Neck Curl Cath Catheter, 2 Cuffs, Right, 62.5 cm     | 1         |
| Kit      | 8888413831 | Swan Neck Curl Cath Catheter Kit, 2 Cuffs, Right, 62.5 cm | 1         |



Swan Neck Oreopoulos-Zellerman Catheters

The Swan Neck Oreopoulos-Zellerman (also called Toronto Western) Catheter is available in Tenckhoff and Missouri configurations. These catheters are differentiated from other Swan Neck catheters by the two silicone discs added to the internal catheter segment.



| Type     | Code       | Description                                                      | Ship Case |
|----------|------------|------------------------------------------------------------------|-----------|
| Catheter | 8888413005 | Swan Neck Oreopoulos-Zellerman Catheter, 2 Cuffs, Left, 44.5 cm  | 1         |
| Catheter | 8888413013 | Swan Neck Oreopoulos-Zellerman Catheter, 2 Cuffs, Right, 44.5 cm | 1         |





## Swan Neck Missouri Catheters

The Swan Neck Missouri Catheter is the same basic design as the Swan Neck Tenckhoff and Curl Cath catheters with the addition of a silicone bead and a felt disk.

| Type     | Code       | Description                                          | Ship Case |
|----------|------------|------------------------------------------------------|-----------|
| Catheter | 8888412601 | Swan Neck Missouri Catheter, 2 Cuffs, Left, 44.5 cm  | 1         |
| Catheter | 8888412619 | Swan Neck Missouri Catheter, 2 Cuffs, Right, 44.5 cm | 1         |




| Type     | Code       | Description                                                   | Ship Case |
|----------|------------|---------------------------------------------------------------|-----------|
| Catheter | 8888413401 | Swan Neck Curl Cath Missouri Catheter, 2 Cuffs, Left, 62.5cm  | 1         |
| Catheter | 8888413419 | Swan Neck Curl Cath Missouri Catheter, 2 Cuffs, Right, 62.5cm | 1         |



## Moncrief-Popovich Swan Neck Curl Cath Catheter


| Type     | Code       | Description                                                      | Ship Case |
|----------|------------|------------------------------------------------------------------|-----------|
| Catheter | 8888414029 | Moncrief-Popovich Swan Neck Curl Cath Catheter, 2 Cuffs, 62.5 cm | 5         |




## Swan Neck Presternal Catheter

**Description of Components:** The presteral (upper) segment of this catheter is tunneled into the subcutaneous tissue of the chest wall. The upper segment connects with the intraperitoneal (lower) segment by a double-barbed titanium connector. Both the presteral and intraperitoneal segments can be trimmed in order to custom fit the patient. The intraperitoneal segment is a variation of the Missouri catheter:


| Type     | Code       | Description                                                                           | Ship Case |
|----------|------------|---------------------------------------------------------------------------------------|-----------|
| Catheter | 8888414011 | Swan Neck Presternal Catheter, 2 Cuffs, 112.8 cm<br>(Swan Neck: 60.3 cm + 52 cm Curl) | 1         |




### Assembled Presternal Catheter



### Presternal Upper Segment 60.3 cm



### Intraperitoneal Lower Segment 52 cm



# Curl Cath Peritoneal Catheters

## Argyle™\* Curl Cath Peritoneal Catheters

### Features

The Curl Cath catheter, a modified Tenckhoff catheter, is suitable for acute or chronic peritoneal dialysis. It can be inserted either percutaneously or surgically.

The Curl Cath catheter is made of translucent, medical-grade silicone rubber with a radiopaque stripe. It has a 4.9 mm O.D., and a 2.6 mm I.D., and comes in adult or paediatric lengths. Unlike polyurethane catheters, the Curl Cath can be extended or repaired if necessary, eliminating the need for costly catheter replacement due to external damage.

Felt cuffs anchor the catheter subcutaneously. To meet the varying needs of patient and physician, two cuff configurations are offered:

- One attached cuff – in the preperitoneal position
- Two attached cuffs – one subcutaneous and one preperitoneal

The Ultem™† adapter included in the packaging provides a tight closure with standard CAPD transfer sets.





Titanium adapters are also available.

† Made with PolyTetra FluoroEthylene (PTFE)

### Product Highlights

- Curled tubing with numerous inflow/outflow holes
- Can be inserted percutaneously or surgically
- Made of translucent, medical-grade silicone rubber
- Radiopaque stripe
- Can be extended or repaired if necessary
- Sterile

### Curl Cath Catheters – One Cuff





| Type     | Code       | Description                    | Ship Case |
|----------|------------|--------------------------------|-----------|
| Catheter | 8811313014 | Curled Catheter, 1 Cuff, 57 cm | 1         |


| Type     | Code       | Description                        | Ship Case |
|----------|------------|------------------------------------|-----------|
| Catheter | 8811313013 | Curled Catheter, 1 Cuff, 60 cm     | 1         |
| Kit      | 8817278010 | Curled Catheter, 1 Cuff Kit, 60 cm | 5         |

| Type     | Code       | Description                                       | Ship case |
|----------|------------|---------------------------------------------------|-----------|
| Catheter | 8888425702 | Oreopoulos-Zellerman Curled Catheter 1 Cuff, 63cm | 1         |


Curl Cath Catheters – Two Cuffs



| Type     | Code       | Description                         | Ship Case |
|----------|------------|-------------------------------------|-----------|
| Catheter | 8811313015 | Curled Catheter, 2 Cuffs, 57 cm     | 1         |
| Kit      | 8817278007 | Curled Catheter, 2 Cuffs Kit, 57 cm | 5         |



| Type     | Code       | Description                         | Ship Case |
|----------|------------|-------------------------------------|-----------|
| Catheter | 8811313010 | Curled Catheter, 2 Cuffs, 62 cm     | 1         |
| Catheter | 8888411702 | Curled Catheter 2 Cuffs, 62 cm      | 1         |
| Kit      | 8817278006 | Curled Catheter, 2 Cuffs Kit, 62 cm | 5         |





| Type     | Code       | Description                            | Ship Case |
|----------|------------|----------------------------------------|-----------|
| Catheter | 8811313009 | Curled Catheter, 2 Mobile Cuffs, 62 cm | 1         |

| Type     | Code       | Description                                                 | Ship Case |
|----------|------------|-------------------------------------------------------------|-----------|
| Catheter | 8888425710 | Oreopoulos-Zellerman Curled Catheter 2 cuffs,1 bubble 62 cm | 1         |


# Paediatric Peritoneal Catheters

## Paediatric Tenckhoff Catheters


| Type     | Code       | Description                         | Ship Case |
|----------|------------|-------------------------------------|-----------|
| Catheter | 8888414219 | Paediatric Tenckhoff, 1 Cuff, 30 cm | 1         |





| Type     | Code       | Description                                    | Ship Case |
|----------|------------|------------------------------------------------|-----------|
| Catheter | 8812329001 | Neonatal Tenckhoff, 1 Subcutaneous Cuff, 31 cm | 1         |




| Type     | Code       | Description                                      | Ship Case |
|----------|------------|--------------------------------------------------|-----------|
| Catheter | 8812321001 | Paediatric Tenckhoff, 1 Subcutaneous Cuff, 37 cm | 1         |




| Type     | Code       | Description                          | Ship Case |
|----------|------------|--------------------------------------|-----------|
| Catheter | 8888414201 | Paediatric Tenckhoff, 2 Cuffs, 30 cm | 1         |




| Type     | Code       | Description                        | Ship Case |
|----------|------------|------------------------------------|-----------|
| Catheter | 8812329002 | Neonatal Tenckhoff, 2 Cuffs, 31 cm | 1         |



| Type     | Code       | Description                          | Ship Case |
|----------|------------|--------------------------------------|-----------|
| Catheter | 8888414227 | Neonatal Tenckhoff, 2 Cuffs, 31 cm   | 1         |
| Catheter | 8810890014 | Paediatric Tenckhoff, 2 Cuffs, 32 cm | 1         |




| Type     | Code       | Description                          | Ship Case |
|----------|------------|--------------------------------------|-----------|
| Catheter | 8810890003 | Paediatric Tenckhoff, 2 Cuffs, 37 cm | 1         |





### Paediatric Curl Cath Catheter

| Type     | Code       | Description                                  | Ship Case |
|----------|------------|----------------------------------------------|-----------|
| Catheter | 8815677001 | Paediatric Curl Cath Catheter, 1 Cuff, 39 cm | 1         |





### Paediatric Tenckhoff Swan Neck Catheter

| Type     | Code       | Description                                                     | Ship Case |
|----------|------------|-----------------------------------------------------------------|-----------|
| Catheter | 8888410506 | Paediatric Swan Neck Tenckhoff Catheter, 2 Cuffs, Left, 37.5 cm | 1         |




### Paediatric Swan Neck Curl Cath Catheters


| Type     | Code       | Description                                           | Ship Case |
|----------|------------|-------------------------------------------------------|-----------|
| Catheter | 8888413100 | Infant Swan Neck Curl Cath Catheter, 2 Cuffs, 38.9 cm | 1         |




| Type     | Code       | Description                                                   | Ship Case |
|----------|------------|---------------------------------------------------------------|-----------|
| Catheter | 8888414813 | Paediatric Swan Neck Curl Cath Catheter, 2 Cuffs, Left, 42 cm | 1         |



| Type     | Code       | Description                                             | Ship Case |
|----------|------------|---------------------------------------------------------|-----------|
| Catheter | 8888413101 | Paediatric Swan Neck Curl Cath Catheter, 2 Cuffs, 43 cm | 1         |





| Type     | Code       | Description                                             | Ship Case |
|----------|------------|---------------------------------------------------------|-----------|
| Catheter | 8888413102 | Adolescent Swan Neck Curl Cath Catheter, 2 Cuffs, 59 cm | 1         |



### Paediatric Swan Neck Curl Cath Missouri Catheter

| Type     | Code       | Description                                                            | Ship Case |
|----------|------------|------------------------------------------------------------------------|-----------|
| Catheter | 8888415307 | Paediatric Swan Neck Curl Cath Missouri Catheter, 2 Cuffs, Left, 42 cm | 1         |



# Accessories For Peritoneal Catheters

## Peritoneal Catheter Straightening Stylet

Keeps the catheter straight for maneuvering during insertion into the peritoneum. The stylet can be used on all catheter types. Non-sterile. Autoclavable.


| Code       | Description                                                                    | Ship Case |
|------------|--------------------------------------------------------------------------------|-----------|
| 8888415661 | Insertion/Straightening Stylet, 61 cm, for use with catheters 61 cm or shorter | 1         |



## Atraumatic Peritoneal Tunneling Stylet

This atraumatic tunneling stylet penetrates tissue cleanly. As the stylet dilates the tissue, the practitioner simultaneously pulls the catheter through the tunnel. Made of pliable stainless steel, it can be bent to create a customized, curved tunnel. Sterile, single use. The stylet is 5 Fr x 23 cm (9").

| Code       | Description                            | Ship Case |
|------------|----------------------------------------|-----------|
| 8817791001 | Atraumatic Peritoneal Tunneling Stylet | 5         |



## Faller Tunneling Stylet

This stylet is an alternative to the atraumatic tunneling stylet. It has a cutting edge on the portion leading into the tunnel and a barbed portion on the opposite end which connects to the catheter. Sterile, single use.


| Code       | Description             | Ship Case |
|------------|-------------------------|-----------|
| 8888415679 | Faller Tunneling Stylet | 1         |



## Titanium Adapter

An alternative to the Beta-Cap Adapter. Can be used on all Argyle™\* peritoneal catheters.

| Code       | Description               | Ship Case |
|------------|---------------------------|-----------|
| 8888415604 | Titanium Two Part Adapter | 1         |



## Argyle™\* Pull-Apart Introducer Set with PTFE

Sterile contents included in the kit:

- |                                      |                                 |
|--------------------------------------|---------------------------------|
| (1) 16 Fr Pull-Apart Sheath/Dilator† | (1) 12 mL Syringe               |
| (1) 18 ga Introducer Needle          | (1) J/Straight 0.038" Guidewire |

| Code       | Description                      | Ship Case |
|------------|----------------------------------|-----------|
| 8815544010 | 16 Fr Pull-Apart Introducer Set† | 5         |

† Made with PolyTetra FluoroEthylene (PTFE).



## Medical Grade Silicone Adhesive Type A

This sterile adhesive is used for:

- Bonding felt cuffs to the catheter
- Repairing damaged catheters
- Extending the length of existing catheters

| Code       | Description                            | Ship Case |
|------------|----------------------------------------|-----------|
| 8810807001 | Medical Grade Silicone Adhesive Type A | 1         |



## Injection Sealing Caps

Used to close off the end of any Argyle™\* Catheter. Universal luer fitting.

| Code       | Description           | Ship Case |
|------------|-----------------------|-----------|
| 8813791001 | Injection Sealing Cap | 50        |



## Beta-Cap Adapter

These adapters form a tight closure with standard CAPD transfer sets, or connect to a sealing cap. The Beta-Cap adapter is made of Ultem™\* resin. Can be used with all Argyle™\* Peritoneal Catheters.

| Code       | Description      | Ship Case |
|------------|------------------|-----------|
| 8814661001 | Beta-Cap Adapter | 25        |



## Beta-Cap Clamp

Used to clamp the peritoneal catheter tubing when not in use.

| Code       | Description    | Ship Case |
|------------|----------------|-----------|
| 8810805001 | Beta-Cap Clamp | 25        |



## Catheter Extender

Provides 15.2 cm of silicone tubing with a titanium connector.

| Code       | Description                      | Ship Case |
|------------|----------------------------------|-----------|
| 8888415612 | Adult Titanium Catheter Extender | 1         |




## Peri-Patch Repair Kit

Use an Argyle™\* Peritoneal Catheter Repair Kit to replace the external portion of a damaged catheter and eliminate the need for catheter replacement. Sterile.

**NOTE:** To repair an indwelling catheter there must be at least 3 cm of undamaged tubing at the exit site.

| Code       | Description           | Ship Case |
|------------|-----------------------|-----------|
| 8888285001 | Peri-Patch Repair Kit | 1         |



**NOTE:** Also order Medical Grade Silicone Adhesive using code # 8810807001.

### Sterile contents included in the kit:

- (1) 15 cm Silicone Rubber Catheter Extension (2.6 mm I.D., 4.9 mm O.D.)  
with Double-Barbed Connector
- (1) Beta-Cap Adapter, Cap and Clamp
- (1) Glue Mold



**Mozarc Medical**  
EMEA Regional Headquarters  
Via Varesina, 162  
20156 Milano MI, Italy

 [emeacontactus@mozarcmedical.com](mailto:emeacontactus@mozarcmedical.com)

 [mozarcmedical.com](https://mozarcmedical.com)

 [LinkedIn](#)

Copyright © 2024 Mozarc Medical Holding LLC.  
Mozarc, Mozarc Medical, the Mozarc Medical logos,  
and Empowering Patients. Enriching Lives.  
are trademarks of Mozarc Medical.

<sup>TM</sup>\* Third-party brands are trademarks of their respective owners.  
\* Argyle<sup>TM</sup> is a trademark of Cardinal Heath, Inc., used under license.

09/2024. EMEA-RC-2300063

Products may not be available in certain countries.

**Mozarc**  
*Empowering patients.  
Enriching lives.* **medical**