

19F Third Lumen
(Mahurkar™* triple
lumen catheters)

Mahurkar™* Elite acute catheters

A full spectrum of acute dialysis care

Laser-etched priming volumes on clamps

Reinforced hub with
catheter sizing label

High flow internal lumen diameter



Double-D™*
lumen design




Modified Double-D™*
lumen design

Visual insertion depth markings

Wrapped, laser-cut
side slots

Symmetrical tip design



Sizing and configuration options

Mahurkar™* 12 Fr acute dual lumen catheter lengths (cm)

High flow: Up to 400 mL/min	Straight Extension	Curved Extension	Pre-Curved Shaft Extension
	13 cm	13 cm	13 cm
	16 cm	16 cm	16 cm
	20 cm	20 cm	20 cm
	24 cm	24 cm	
	30 cm		

Mahurkar™* 13 Fr acute dual lumen catheter lengths (cm)

High flow: Up to 450 mL/min	Straight Extension	Curved Extension	Pre-Curved Shaft Extension
	13 cm	13 cm	13 cm
	16 cm	16 cm	16 cm
	20 cm	20 cm	20 cm
	24 cm	24 cm	
	30 cm		

Mahurkar™* 12.5 Fr acute triple lumen catheter lengths (cm)

Third lumen 19G High flow: Up to 400 mL/min	Straight Extension	Curved Extension
	13 cm	13 cm
	16 cm	16 cm
	20 cm	20 cm
	24 cm	24 cm
	30 cm	

Learn more: View full product details in the Mahurkar™* Elite acute catheter brochure.

©2024 Mozarc Medical Holding LLC. Mozarc, Mozarc Medical, the Mozarc Medical logos, and Empowering Patients. Enriching Lives. are trademarks of Mozarc Medical. Covidien and Medtronic are trademarks of Medtronic used under license. MAHURKAR™*, Double-D™*, and the Double-D design are U.S. registered trademarks of Sakharam D. Mahurkar, used under license.

™* Third-party brands are trademarks of their respective owners. 02/2024 - US-RC-2400042 – [IH]