

Ordering Information Chart

Palindrome™ Precision chronic catheters

Product Codes Kits	Product Codes Sport Packs	Insertion Length	Overall Length
Palindrome™ Precision catheter			
8888145014P	8888119360P	19 cm	36 cm
8888145015P	8888123400P	23 cm	40 cm
8888145016P	8888128450P	28 cm	45 cm
8888145017P	8888133500P	33 cm	50 cm
8888146044P	n/a	44 cm	61 cm
8888145018P	n/a	55 cm	72 cm
Palindrome™ Precision catheter with VenaTrac™ insertion stylets			
8888145039P	8888119364P	19 cm	36 cm
8888145040P	8888123404P	23 cm	40 cm
8888145041P	8888128454P	28 cm	45 cm
8888145042P	8888133504P	33 cm	50 cm
Pre-Curved Palindrome™ Precision catheter			
8888145058P	8888119370P	19 cm	36 cm
8888145059P	8888123410P	23 cm	40 cm
8888145060P	8888128460P	28 cm	45 cm
8888145061P	8888133510P	33 cm	50 cm
Palindrome™ Precision H chronic catheter - Heparin coating			
8888145043P	8888119365P	19 cm	36 cm
8888145044P	8888123405P	23 cm	40 cm
8888145045P	8888128455P	28 cm	45 cm
8888145046P	8888133505P	33 cm	50 cm
Pre-Curved Palindrome™ Precision H chronic catheter - Heparin coating			
8888145068P	8888119371P	19 cm	36 cm
8888145069P	8888123411P	23 cm	40 cm
8888145070P	8888128461P	28 cm	45 cm
8888145071P	8888133511P	33 cm	50 cm

Product Codes Kits	Product Codes Sport Packs	Insertion Length	Overall Length
Precision SI chronic catheter - Silver ion sleeve			
8888145062P	8888119368P	19 cm	36 cm
8888145063P	8888123408P	23 cm	40 cm
8888145064P	8888128458P	28 cm	45 cm
8888145065P	8888133508P	33 cm	50 cm
8888146144P	n/a	44 cm	61 cm
8888145066P	n/a	55 cm	72 cm
Palindrome™ Precision HSI chronic catheter - Heparin coating and silver ion sleeve			
8888145057P	8888119369P	19 cm	36 cm
8888145048P	8888123409P	23 cm	40 cm
8888145049P	8888128459P	28 cm	45 cm
8888145050P	8888133509P	33 cm	50 cm

More product details

View the Palindrome™ Precision chronic catheter brochure

[Click](#) or scan to view



Components

Components	Qty	Sport Packs	Kits
Palindrome™ Precision base, H, SI, and HSI catheters			
14.5Fr symmetric tip catheter	1	•	•
VenaTrac™ over-the-wire insertion stylets	2	•	•
16Fr Valved pull-apart safety sheath	1	•	•
Bifurcated tunneler	1	•	•
12Fr Tissue dilator	1	•	•
14Fr Tissue dilator	1	•	•
Injection sealing caps	2	•	•
18g Introducer needle	1		•
J/Straight 0.038" guidewire	1		•
12mL Syringe	1		•
#11 Scalpel	1		•
Telfa™ island dressings	2		•
4" x 4" cotton gauze sponges	4		•

Repair Kits

Repair Kits	Product Codes	Qty
Hemodialysis Catheter Repair Kit	8888200001	1

Components

Hemodialysis Catheter Repair Kit	Qty
Arterial repair assembly	1
Venous repair assembly	1
Sealing caps	2
Temporary slide clamps	2
Drape	1
Disposable scissors	1
Measurement guide	1
Priming volume label sheet	1

Color Key

- Palindrome™ Precision chronic catheter (base)
- Palindrome™ Precision H chronic catheter (heparin coated)
- Palindrome™ Precision SI chronic catheter (silver ion sleeve)
- Palindrome™ Precision HSI chronic catheter (heparin coated and silver ion sleeve)

710 Medtronic Parkway
Minneapolis, MN 55432
800.962.9888
508.261.8000
MozarcMedical.com